



# Bodystation II

**High Performance**

Photo System for Aesthetics

# Aesthetics **Bodystation II**

## High Performance **Photo System** for Aesthetics

### **The Mobile Photo Studio for Consistent Results**

The Bodystation II is our most advanced mobile imaging solution for consistent before and after photo documentation of virtually any aesthetic face or body procedure. The redesigned Bodystation II can operate “unplugged” for several hours\* for even greater mobility within your office, eliminating the need for a dedicated photo room.

### **Show the Success of Your Treatments**

Patients often accommodate to the outcome of their treatment results over time. The Bodystation II reinforces the quality of aesthetic results by accurately tracking the progress of procedures. Stored images are easily accessible and can be retrieved by the FotoFinder query feature. Whether you are marketing to prospective clients or tracking the effectiveness of procedures, you can showcase side-by-side photo comparisons and print impressive reports.

### **Ease Of Use is Key**


Computer controlled camera settings and clear onscreen instructions make it easy to take excellent photos. FotoFinder’s Laser Liner displays a red laser line for repeatable patient positioning and accurate photo capturing eliminating the need for a floor mat. No photography experience is necessary for system use and only little amount of staff training is required.



\*Battery runtime may vary depending on selected laptop power plan. System must be plugged in for printing.



## Aesthetics **Bodystation II**




### FotoFinder Guided Photography Templates

assist to quickly capture full face or body images. Choose from a large variety of photography templates that display step-by-step onscreen instructions. FotoFinder's virtual patient shows your staff treatment-specific patient poses and body area to be photographed. Every photo is automatically labeled and stored in your database for effortless retrieval.



### FotoFinder Ghost

feature displays a transparent baseline photo on top of the live image, making it easy to align your patient at follow ups. Onscreen instructions, consistent camera settings, Laser Liner and FotoFinder Ghost make follow ups from the right angle and same distance a cinch.



 **FotoFinder**<sup>®</sup>



### FotoFinder Zoom Ring

When taking photos of face or body, the focal length (zoom) on the camera lens must be adjusted. The Zoom Ring makes it easy to find the right setting with its precision integrated stopping points and color coding. FotoFinder software verifies the correct setting is being used, providing immediate quality control.



### FotoFinder Proposer

is an impressive client consultation tool for effective, time-saving customization. Simply highlight treatment areas and drag and drop your recommended procedures on your client's photo. Add from your personalized product database your recommended skin care products and treatment protocols and present your proposal on a professional report for an exciting patient experience.

### FotoFinder UVscan

module converts digital color photos into impressive monochrome images helping providers to promote their rejuvenation treatments effectively. It's simple, fast and does not require additional UV-photo equipment.



## FotoFinder Universe iPad App

offers access to your images from anywhere in your office, for viewing, editing and comparison. Quick snapshots can be taken and are automatically archived in your FotoFinder database. Create treatment plans on your iPad during consultation and review the results with your clients to reinforce treatment success.

## Compelling Marketing and Reports

FotoFinder's time saving editing tools allow you to quickly create impressive side-by-side comparisons for your online marketing activities. Choose from a variety of compelling print layouts, which feature your company logo and address for a polished and professional look. Improve your social media presence and get more referrals of satisfied patients.

## Aesthetics Bodystation II

### High Performance Photo System

#### Elegant Ergonomic Design

The Bodystation II looks elegant and features a convenient drawer for storage of printing paper and an extractable mouse pad for left and right-handed use.

#### Laptop

Preinstalled with FotoFinder Universe software for a hassle-free out of box experience.

#### Battery Pack

FotoFinder's built-in battery allows several hours of mobile image capturing for greater mobility\*.

#### Large Lockable Castors

4" lockable castors roll easily over tile, laminate, hardwood, or carpet.

#### High resolution DSLR Camera

24 MP camera tethered to the laptop and controlled by FotoFinder Universe software. No importing from an SD card required.

#### Internally Wired Cables

No loose cables for tangle-free efficient use.

#### Laser Liner

Built-in laser projecting a red line for patient positioning eliminating old fashioned floor mats.

#### Photo Printer

Create a photo album for your waiting room and give your client a professional consultation or before & after report.

\*Battery runtime may vary depending on selected laptop power plan. System must be plugged in for printing.

Contact us to schedule a live online demo!

**888-501-0805**

E-mail: [info@fotofinder-systems.com](mailto:info@fotofinder-systems.com)

Phone: +1-443-283-3865

**[www.fotofinder-systems.com](http://www.fotofinder-systems.com)**