



Bodystudio

# ATBM

master

4<sup>th</sup> Gen

Automated Total  
Body Mapping



[www.fotofinder-systems.com](http://www.fotofinder-systems.com)

©FotoFinder Systems GmbH  
USA - V1.0 / 2024-03

# ATBM master

4<sup>th</sup> Gen

The **ATBM master 4th Generation** is the latest and most refined skin imaging system FotoFinder developed ever. Building on years of experience we finetuned its technology and pushed the limits further to build the most sophisticated and reliable FotoFinder machine for Total Body Dermoscopy.

FotoFinder's special RAW image processing and cross polarized camera system with laser focus allow skin imaging at spectacular image quality. Experience a new, streamlined, and intuitive interface and bring Total Body Mapping to the next level.

## Intelligent, Multifunctional Software Platform Universe 3.5

Versatile in its use, with modules for medical imaging, dermoscopy, trichoscopy and aesthetics.

## High-resolution Imaging with PolFlash XE

Camera system with proprietary xenon flash. Whole body photos in RAW format, cross-polarized or non-polarized. Photo and video recording for aesthetics.

## Silent Medical Server with Powerful Q-Processor

The Q-Processor, the core of the ATBM master, ensures exceptional performance combined with massive storage capacity and built in data redundancy.

MDR approved and licensed by Health Canada as Class II device.



## True Video Dermoscopy with medicam 1000s

Seamlessly integrated Full HD video camera featuring CrystalView technology and optical live zoom.

## Faster Total Body Dermoscopy Workflow

Streamlined workflow for capturing and reviewing images significantly reduce examination time.

## Aesthetics Module with Guided Photography and Ghosting

Consistent images and videos of face and body treatments, Proposer for treatment planning and skin care recommendations, visualization of redness and sun damage and skin texture – all in one platform included.

## Turnkey Installation

Pre-installed ATBM tower with automatic camera lift and FotoFinder's high-performance hardware, suitable for single or multi-user installations.

## Mosaic View

### Unmask the “Ugly Duckling”

The **Mosaic View** displays detected lesions from ALL total body images and intelligently organizes them on a single screen. This Mosaic view can help providers identifying outlying lesions.

During image comparison the **Bodyscan** can assist the provider identifying changes in existing lesions and new ones through differently colored circles.





# MEDICAM 1000s

## medicam 1000s V4

A Benchmark for Video Dermoscopy Worldwide.

### Performance and quality redefined:

The reengineered **medicam 1000s** delivers sharper images and enhanced optical zoom, with refined software control for even faster performance. The ultimate allrounder for skin, hair and nails.



### Brilliantly Engineered

Crafted from high quality, anodized, lightweight aluminum, the black medicam body offers both durability and sophistication.

### Brilliant and Reliable

Designed for daily use and speedy operation, the camera and lens are extremely durable yet flexible, ensuring an ergonomic workflow.

### Precision Optics

The **D-Scope IV** lens with TwinLight technology allows easy toggling between polar-ized and non-polarized lighting..

### Excellent Clinical Images

The floodlight illumination with computer-controlled LEDs illuminate skin, scalp, and nails for outstanding clinical pictures. For consistent macro imaging optional magnetic spacers are included.

## D-Scope IV

A New Vision in Dermoscopy.



The **D-Scope IV** precision lens stands for exceptional image quality. Equipped with optics boasting high luminous intensity, **D-Scope IV** provides crystal-clear dermoscopy images with variable optical zoom and autofocus. Its impressive live zoom allows real-time examination with unrivaled detail.



Polarized Contact Dermoscopy



Polarized Non-contact Dermoscopy



Clinical Image



## Features at a glance

### Spectacularly Sharp

- ▶ Full HD CrystalView technology
- ▶ Brilliance in dermoscopy of skin, hair and nails

### Intelligent TwinLight Technology

- ▶ Toggle between polarized and non-polarized dermoscopy
- ▶ Convenient toggle button on the medicam control panel

### Precision Lens System

- ▶ Optical Live-zoom
- ▶ Variable magnification from 20× to 140×
- ▶ Autofocus on every zoom level

### Contact or Non-contact Dermoscopy

- ▶ Quick change front tips
- ▶ Dermoscopy with and without glass contact plate
- ▶ Conical tip for lesions between fingers and toes



Made in Germany. With love.



# Master in Aesthetics

## Impressive Consultations and Before and After Images.

With the **ATBM master** you have a versatile, high-end imaging platform for consistent imaging of face and body. Guided Photography templates, automated optical zoom, and FotoFinder's Ghost feature for follow ups are key elements for consistent before and after imaging. Present your artistry on social media with spectacular side by sides slides or concurrent comparison video clips created on your machine.

**Proposer V10** is your consultation platform to visualize pigmentation, redness, or skin texture and to create convincing treatment and skin care reports. Boost your consultations and convert more clients with impressive, branded Proposer reports.



# FotoFinder Platinum Care



## We keep your FotoFinder in shape.

**FotoFinder Platinum Care** ensures optimal performance of your **FotoFinder ATBM Master**. As an investment critical to the success of your practice, you receive priority support from our friendly US-based team. We promptly address any issues or questions your staff may encounter. Our knowledgeable technicians install software updates at your convenience, and in case of hardware problems, loaner items are dispatched promptly to minimize downtime. Rely on FotoFinder's dedicated service team for personalized support, including training and assistance with operations or complex network installations.

# Hospital System Integration

## Client server installation and EMR integration.

FotoFinder **Universe** software is a flexible client server application which can be easily deployed in large clinic networks.

Multiple camera systems operate simultaneously with one or more databases. Image access through virtual clients enables access to FotoFinder data anywhere in your clinic. Optional integration through HL7 and DICOM is available to transmit patient information and images automatically between FotoFinder and the EMR.





**FotoFinder Systems, Inc.**

7100 Columbia Gateway Dr, Suite 120  
Columbia, MD 21046, USA

Toll-free: +1 888 501 0805  
Phone: +1 443 283 3865  
Fax: +1 443 283 3869  
Mail: [info@fotofinder-systems.com](mailto:info@fotofinder-systems.com)

For more information visit

**[www.fotofinder-systems.com](http://www.fotofinder-systems.com)**